

White Paper

The Role of PCR in Molecular Biology

By Cedric Herry Director of Research and Development, Erlab Inc.



In 2020, while the world suffered from the worse pandemic in decades, the acronym PCR, became familiar to most of us. PCR (Polymerase Chain Reaction) procedure has recently been realized as a major tool specifically in diagnosing people contaminated by Coronavirus. PCR was invented in the early 80's and became one of the most widely used methods in molecular biology. It is used mainly in labs to produce sufficient quantities of DNA or RNA from a sample where this material is in limited quantities.



The ultimate goal of PCR is to produce multiple identical DNA fragments from one single DNA or RNA strand. PCR is often used in multiple applications such as forensics, detection of microorganisms, prenatal diagnosis, medical diagnostics, and agriculture, to name a few. See the chart to the right for an example of reverse transcription polymerase chain reaction which is a laboratory technique combining reverse transcription of RNA into DNA, along with amplification of specific DNA targets using PCR.



This contamination can come from the environment, or from cross contamination by other preamplified samples. Depending on the application, the impact of contaminations can be dramatic and could result in bad diagnostics or could compromise criminal investigations.

In order to optimize PCR and to avoid compromised results, the PCR process must be clearly organized, following 3 major steps, each of them requiring stringent protocols:

- 1 Lysis of cells and extraction of DNA or RNA
- 2 Preparation of PCR mixture before amplification
- 3 Amplification by itself in a thermocycler

DNA and RNA are contained inside cells, and more generally, in the cell nucleus. Prior to PCR, the first step is the extraction of DNA or RNA from the cell. This can be done using different strategies. During this step, the sample must be protected to avoid it being contaminated by environmental pollution, while at the same time the operator must be protected against hazardous chemical that are used to extract DNA or RNA. If the sample is a pathogenic micro-organism, this part of the procedure must be performed in a Biological Safety Cabinet type 2. At the end of step one, the extracted sample must be isolated, and the next step must absolutely be performed in a second work station.

During step 2, also called pre-PCR, extracted samples are mixed with all reactants that are needed to perform amplification. These reactants include primers that hybridize into flanking sequences on opposing strands of the target, 4 deoxyribonucleoside triphosphates and a DNA polymerase along with buffer, co-factors of enzymes and water. During this sequence, the sample is highly sensible to both environmental and cross contaminations, as DNA or RNA strands are no longer protected by cells, and are in very limited quantities. Environmental contamination could come from dust air pollution that embeds micro-organisms or other biological material, and cross-contamination could come from a previous sample that was manipulated in the same enclosure.



The best way to protect the sample against air contamination, is to blow purified air on the sample. Air is purified by passing through high efficiency particulate filters such as HEPA H14 or ULPA U16. These filters are able to trap airborne pollution due to a network of non-woven glass fibers, that traps submicronic particles, at very high efficiency (more than 99.995% for H14 filters according to EN 1822-1 or 99.9995% for ULPA U16 filters according to the EN1822-1). When such filters are correctly used, and sized appropriately, ISO 5 air cleanliness is accomplished.



PCR Role in Molecular Biology



Achieving such air quality, requires the use of appropriate filters, to correctly manage airtightness through the filter, and to deliver enough air flow to maintain a sufficient overpressure when compared to the surrounding environment. Working in this type of purified air is mandatory for PCR applications. It can be achieved by working in a cleanroom, but more commonly it is done through the use of dust-free enclosures, using HEPA or ULPA filters.

Additionally, during step 2, it is necessary to prevent cross contamination. The risk is high to contaminate a fresh extracted sample with a sample previously mixed in the enclosure. To avoid such contamination, UV-C, at 254 nm has been shown to be efficient. Short UV waves at 254 nm directly affect cellular RNA and DNA (e.g. bacteria and viruses) and degrade them. UV absorption creates bonds between adjacent nucleic acids. These alterations in DNA or RNA block the possibility of replication. When buying a safety enclosure for pre-PCR, it is necessary to ensure that the proper wavelength is delivered by a UV bulb and that it is delivering a UV dose strong enough to destroy contamination between the 2 preparations. For that, the ratio between the workstation dimension, the UV power, and its distance to the worktop are some of the

most important parameters. As an example, RNA viruses are deactivated by 12 mJ/cm2 of UV-C lightning according to Predicted Inactivation of Viruses of Relevance to Biodefense by Solar Radiation, *C. David Lytle and Jose-Luis Sagripanti**, *J Virol. 2005 Nov*; 79(22): 14244-14252.

If these first 2 steps are correctly performed, then step 3, PCR by itself can be performed in a thermocycler. As the samples are sealed, the risk of contamination in thermocyclers is close to zero.

For more than 2 decades, Erlab has delivered thousands of PCR workstations. The most recent models are dedicated to 2 types of situations:

CaptairBio 320 is recommended when pre-PCR is performed in a cleanroom. The CaptairBio 320 Smart is a workstation where UV contamination is optimized to avoid cross-contamination. It features an 18W UV bulb at 254 nm that delivers a minimum of 0.08 jJ/s/cm2 at the level of the worktop.

CaptairBio 321 and 391, are 2 workstations that prevent both cross-contamination and environmental contamination of samples. They can both embed HEPA or ULPA filters and are equipped with high performance 254 nm LED 900/950 lux UV bulbs.





About the author

Dr. Cedric Herry is Erlab's Director of R&D, a position he has held for many years. He is an expert in air filtration and safety and sits on the European Standardization Committees for laboratory equipment. Dr. Herry's specialities include: Air treatment, Chemical Hazards, and Sustainable Laboratories.



The Erlab Research and Development laboratory

About Erlab

Since 1968, Erlab has been a specialist, inventor and world leader in ductless, zero-emission filtering fume hoods for laboratories to provide total safety in chemical handling.

Erlab filtration

We provide technologies to protect laboratory staff from inhaling chemicals. This is made possible thanks to our Research and Development (R&D) department, which has continuously improved our filtration technology for more than 50 years. That's why, in 2009, we invented the ERLAB ABOVE label for tried and tested filtration technology.

The AFNOR NF X 15-211: 2009 standard

Erlab's filtration technology conforms to the NF X 15-211: 2009 standard, the industry's most demanding standard for molecular filtration, developed by a committee of independent scientists and specialized manufacturers.

This text imposes performance criteria linked to:

- Filtration efficiency
- Containment efficiency
- Air face velocity
- · Documentation: chemical listing

The ESP program

A set of three services included with the purchase of each device designed to ensure your safety.

eValiQuest Risk analysis – Determination of protection needs – Determination of ergonomic needs.

- ValiPass
- Certified installation Total safety for handling.



Ongoing monitoring - Preventative and maintenance inspections - Device reconfiguration based on protection needs - Development of handling.

Flex technology

The combination of molecular and particulate filtration technologies allows a single device to meet laboratories' protection needs. This innovation from Erlab's R&D department offers unprecedented flexibility, versatility and value. A single device can be reconfigured over time and easily reassigned to other applications.

Smart technology

Smart technology is a simple and innovative means of communication that improves safety. This technology uses a light and sound signal to indicate the user's level of protection. The advantages of the technology are:

1/ Light pulsation: Real-time communication via LED light pulses intuitively alerts the user to the device's operating status.

2/ Simplicity: One-touch activation.

3/ Detection system: The exclusive detection system continuously monitors filtration performance.

4/ Built-in monitoring: This service provides direct access to the status, settings and history of your device.

France

Spain

United Kingdom +44 (0) 1722 341 940 | export.north@erlab.net +39 (0) 2 89 00 771 | export.south@erlab.net

